

Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

FAST SHIPPING AND DELIVERY

- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED
 SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

WE BUY USED EQUIPMENT


Sell your excess, underutilized, and idle used equipment We also offer credit for buy-backs and trade-ins **www.artisantg.com/WeBuyEquipment** *>*

Instra View REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at **www.instraview.com** *>*

LOOKING FOR MORE INFORMATION?

Visit us on the web at <u>www.artisantg.com</u> for more information on price quotations, drivers, technical specifications, manuals, and documentation



Operating Instructions Mixer Mill Type MM 301



184VA/0,5A; 30%ED; 2x Fuse T2A

Power 120V 50/60Hz

© Retsch GmbH, 42781 Haan, Rheinische Str.36, Germany, Doc.no. D 98.741.9999 Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Notes on the operating instructions2				
Safety4Notes on safety4Warning notes5Repairs5Confirmation5				
Technical Data6Use for the intended purpose6Safety systems6Description of function7Drive7Vibration frequency of the grinding jar fastening7Rated power7Emissions7IP rating7Dimensions7Installation area required7				
Transport and installation8Packaging8Transport8Temporary storage8Scope of supply8Parameters for the installation point8Installation9Electrical connection9				
Operation 10Controls and operation10Graphical display of the controls:10 Table of controls 11Clamping the grinding jar by means of the ratchet fastening.12Grinding jar types12Switching on / off13Setting the time13Setting the frequency13Programming the MM 30114Calling up a program14Adjustment lockout15The grinding process15Changing fuses15				
Notes on operation				
General16The grinding process16Process material quantities and particle sizes16Cleaning17Maintenance17Necessary safety tests17Copyright17Modifications17				
Accessories 18 Further accessories on request. 18 Troubleshooting list 18 Safety regulations (table) 19 Warranty conditions 20				

Warning notes

We use the following signs to warn you against:



Personal injury

Damage to property

Repairs

These operating instructions do not contain any instructions on repairs. For your own safety, only have repairs performed by Retsch GmbH or by an authorised representative (service technicians). 200

In such cases, please notify:

The Retsch agency in your country

Your supplier

Retsch GmbH direct

Your service address:

Confirmation

I have read and noted the chapters entitled Foreword and Safety

Operator's signature

Service technician's signature

Description of function

The milling process in the MM 301 laboratory mill is performed by the beating effect of the grinding balls on the process material and the friction between the grinding balls and the walls of the grinding jars.

The grinding jars, positioned horizontally, are subjected to radial oscillations which ensure that the grinding balls impact against the material one after the other.

The speed can be set from 3 - 30 Hz (120 - 1800 oscillations/minute) and is kept constant by a speed controller during operation.

The milling and mixing duration can be preset digitally from 10 sec. to 99 min.; the setting is preserved in standby mode for further operations.

54

Drive

Single phase AC motor

Vibration frequency of the grinding jar fastening

3 - 30 1/s = 180 - 1800 1/min

The maximum frequency cannot be reached with every loading status: e.g. for a 50 ml steel grinding jar with one grinding ball of 30 mm diameter, only approx. 25 1/s or 1500 1/min can be reached. This overload is indicated by the display alternating between the set and the actually frequency reached.

The MM 301 always operates with the lowest displayed value.

Rated power

Approx. 184 VA

Emissions

Noise measurement to DIN 45635-031-01-KL3 The noise levels are also influenced by the properties of the material to be milled. Example 1: Sound power level $L_{WA} = 71.4 \text{ dB}(A)$ Workplace related emission level $L_{pAeq} = 61 \text{ dB}(A)$ Operating conditions : = 2 steel grinding jars, 25 ml Vessels Milling unit = 1 steel ball, 20 mm each Milled material = crushed quartz approx. 4.0 - 6.0 mmMilled quantity = 8 mlExample 2: Sound power level $L_{WA} = 76 \text{ dB}(A)$ Workplace related emission level $L_{pAeq} = 65 \text{ dB}(A)$ Operating conditions: Vessels = 2 steel grinding jars 5 ml Milling unit = 2 tungstan carbide balls, 8 mm each Milled material = crushed quartz approx. 1.0 - 2.0 mm Milled quantity = 1.5 ml

IP rating

IP40 or IP20 casing aperture for grinding jar fastening

Dimensions

Height: up to approx. 225 mm, width: 300 mm, depth: 470 mm with cover open Height: up to approx. 410 mm, width: 300 mm, depth: 470 mm

Weight: approx. 20 kg, without grinding jars

Installation area required

300 mm x 470 mm; no safety clearances are necessary.

15.03.2006

C Retsch GmbH

7

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Humidity:

Maximum relative humidity 80% at temperatures up to 31 °C, declining in a straight line to 50% relative humidity at 40 °C

Hig mec unk

Higher humidity can cause damage to the electrical and mechanical components, and performance data can change to an unknown extent.

<u>Site altitude:</u>

Max. 2000 m above sea level


Installation

Always set the MM 301 on a firm base.

Remove the braces fitted for transport:

Two such braces **P** are located on the bottom of the MM 301. **Fig.1** They are marked by arrows.

- Loosen bolts S2, see Fig.2
- Loosen the bolts **S2** of the transport brace **P** with washer **S1** by shifting the slot and removing it
- Keep the braces for use if the machine will have to be transported again.



Р

P 52 51

Fig.2

 ∇

If the machine is operated with the transport braces still in place or transported with the braces not in place, mechanical components may be damaged.

The machine must be positioned fully on the laboratory bench. Its cover must not project over the front of the bench.

Otherwise, it may be possible to reach into the machine.

Danger of injury to the hands!

Electrical connection

- See the type plate for the voltage and frequency of the MM 301.
- Ensure that the values agree with those of the available power supply.
- Connect the MM 301 to the power supply using the power cable supplied.

Failure to observe the values on the type plate can lead to damage to electrical and mechanical components.

		1 401	e of controls
Item	Element	Illustration	Function
Α	main switch		Connects the MM 301 to, and isolates it from, the mains
A1	fuse tray		
		í	
В	Power socket	<u>oo</u>	Accommodates the mains cable
С	Ratchet fastening		Accommodates the corresponding grinding jar
		ans a	
D	Hand wheel with		Clamps the grinding jar
D	compression disk		Champs the grinting far
	compression usk		
\mathbf{E}	Locking bolt		Secures against unintended loosening of the hand wheel
F	START button	0	Starts the grinding process
		START	LED lights up
			Access to adjustment mode
G	STOP button		Stops the grinding process before the set time has expired.
		0	LED lights up
		STOP	Press the STOP button a second time = standby
			Access to pause function
н	Display of machine	min	Indicates the pre-selected running time from
	running time	PPP	10 sec. to 99 min
		time	
		time	
J	Rotary knob		Facilitates programming of the desired running time.
	for time setting		· · · · ·
L	Frequency display	1/s	Displays the pre-selected frequency and the operating
		888	frequency.
		frequency	
K	Rotary knob for frequency	A	Regulates motor speed and accordingly the vibration intensit
	display	\cap	of the grinding jars.
	1 2		
Μ	Memory function		Facilitates storage and call-up of up to three times and
			frequency combinations.
	<u> </u>	memory	
Ν	Adjustment lockout	0	Prevents adjustment of time and frequency.
		Č	Press button N to activate the adjustment lockout: the LE
		(A)	will then light up.
			To deactivate the lockout, press button N in standby mode
0	Plexiglas cover	Not shown	Prevents inadvertent reaching into a running machine
P	Transport brace	Not shown	Protects the internal componets of the machine during
L			

15.03.2006

C Retsch GmbH

Doc.no. D 98.741.9999

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com




Fig.2




Fig.3

Switching on / off

The main switch A is located at the rear of the MM 301 Fig.2

- Turn on the main switch A.
- The LED above STOP F flashes briefly after the machine has been turned on. Fig.3
- The time display shows the time last set.
- The frequency display shows the speed last set.

The MM 301 is now in programming mode.

Setting the time

When the MM 301 is switched on it automatically goes into programming mode. The running time last used is shown on display **H**.

Set time to 10 sec - 99 minutes: Fig.3

- By turning knob J counterclockwise you can reduce the time set down to 10 seconds.
- By turning knob J clockwise you can increase the time set up to 99 minutes.

Running time can be set freely between 10 sec. and 99 min. From 10 min. onwards, the time can only be set in full minutes. Turn knob slowly to move in steps of 1 sec. / 1 min.

Turn knob rapidly to move in steps of 30 sec. / 5 min.

Setting the frequency

When the MM 301 is switched on it automatically goes into programming mode. The frequency last used is shown on display L.

Set frequency to 3 - 30 Hz: Fig.3

- By turning knob K counterclockwise you can reduce the frequency set down to 3 Hz.
- By turning knob K clockwise you can increase the frequency set up to 30 Hz.

The frequency can be set freely between 3 and 30 Hz [1/s].

Turn knob slowly to move in steps of 0.1 Hz [1/s]. Turn knob rapidly to move in steps of 5 Hz [1/s].



The MM 301 is equipped with an electronic overload indicator. If there is an overload, the frequency display changes back and forth between the frequency that you have set and the maximum frequency reached.

Lower the frequency to the maximum frequency reached!






Fig.6





The MM 301 has an adjustment lockout N. Button N can be pressed when programming or during operation of the machine. Fig.6

- Press the adjustment lockout N.
- The LED lights up.

It is now impossible to change the time or frequency settings and no values saved to memory can be called up.

The lockout can only be deactivated by pressing button \mathbf{N} in standby mode.

The grinding process

Clamp your grinding jars filled with the process material and the grinding ball/s into the grinding jar fastening. Fig.7

- Close the Plexiglas cover.
- Set the grinding parameters.

Start the MM 301.

Do not open the Plexiglas cover during the grinding process.

Although the grinding jars are brought to a standstill immediately by the builtin brake, the grinding process cannot be continued for the remaining time.

The MM 301 has to be restarted. The initial parameters are then once again available.





Fig.7



Changing fuses

 Two glass fuses (5x20 mm) are required.

 2 fuses mT 0.63 A
 220-240V/50-60Hz

 2 fuses F 1.6 A
 100 / 110 and 120V 50-60Hz

Changing fuses: Fig.8

- Disconnect the mains plug
- Pull out fuse tray **B**
- Replace the fuses
- Slide the fuse tray back in

Cleaning



Do not clean the MM 301 under running water. Danger of fatality from electric shock. Use a moist cloth only. Use no organic solvents.

The ingress protection rating of the MM 301 is IP40 / IP20 (casing aperture for grinding jar fastening).

Maintenance

The MM 301 is to the greatest possible extent maintenance-free. When used in accordance with the manufacturer's instructions the machine requires no maintenance or adjustment.

Necessary safety tests

The limit switch and the brake must be checked for correct functioning every six months.

1. Limit switch

- Remove the analysis sets or grinding jars from their fastenings.
- Switch the MM 301 vibration mill on.
- Set the frequency to 30Hz.
- Start the MM 301.
- Open the Plexiglas cover.
- When the Plexiglas cover has been opened by max. 8 mm, measured from the bottom edge of the housing, the motor brake must be activated.
- The motor must not restart when the Plexiglas cover is closed.

2. Brake

- Switch on the MM 301 vibration mill.
- Set the frequency to 30Hz.
- Switch the MM 301 vibration mill on.
- Open the Plexiglas cover.
- The brake must act within 1 second.

If the machine fails to behave as described in points 1 and 2 above, contact our service technicians.

Copyright

Dissemination or duplication of this documentation, utilisation and dissemination of its contents will only be permitted with the express approval of Retsch GmbH. Contravention will give rise to liability for damages.

Modifications

The right is reserved to modify any specification without notice.

Safety regulations (table)	for the MM 301 from the	for the MM 301 from the individual chapters		
Subject	Action	Hazard		
Safety	Failure to observe the notes on safety can lead to damage to property and personal injury.	We accept no liability whatsoever for any damage or injury in this case.		
Use for the intended purpose	The sample quantity should not fall below 25 % of the grinding jar volume	Otherwise the grinding balls could damage the grinding jar.		
	The materials used for the grinding jar and grinding balls must always be identical.	Otherwise the grinding balls could damage the grinding jar.		
	Do not modify the machine in any way, and use only spare parts and accessories approved by Retsch.	Otherwise the conformity with European Directives declared by Retsch will be invalidated. Furthermore, this will lead to all warranty claims		
	Do not clamp ceramic grinding jars for the	being rendered null and void. If clamped in the grinding jar holder of the MM		
· · · ·	MM 2/200 or jars without complete steel jacket.	301, the metallic parts can cause damage to the ceramic parts.		
Packaging	Please keep the packaging for the duration of the guarantee period.	Your guarantee claim may be invalidated if there is a complaint and the unit is returned in inadequate packaging.		
Transport	The MM 301 must not be jarred, shaken or thrown during transport.	This could damage the electronic and mechanical components.		
Temperature fluctuations	When subjected to major temperature fluctuations, the MM 301 is to be protected from condensation.	Damage to the electronic components may otherwise result.		
Scope of supply	In cases of incomplete delivery and/or transport damage, you must notify the forwarding agent and Retsch GmbH without delay (within 24 hours).	It may not be possible to accept later complaints.		
Ambient temperature	Not below 5 °C. Not above 40 °C.	Higher humidity can cause damage to the electrical and mechanical components, and performance data can change to an unknown extent.		
Humidity	Do not exceed a maximum relative humidity of 80% at temperatures up to 31 °C, declining in a straight line to 50% relative humidity at 40 °C.	Higher humidity can cause damage to the electrical and mechanical components, and performance data can change to an unknown extent.		
Transport braces	Remove the transport braces before commissioning. Secure the machine with transport braces when it is to be transported.	If the machine is operated with the transport braces still in place or transported with the braces not in place, mechanical components may be damaged.		
Setting up the MM 301	The machine must be positioned fully on the laboratory bench. Its cover must not project over the front of the bench.	Otherwise, it many be possible to reach into the machine. Danger of injury to the hands!		
Electrical connection	Ensure that the values on the type plate agree with those of the available power supply.	Electronic components could otherwise be damaged		
Clamping the grinding jar	Do not load the MM 301 unequally.	In order to ensure smooth running of the MM 301, both milling points must be loaded with the same mass.		
	Do not clamp grinding jars from the MM 2/200 or without complete steel jacket.	If clamped in the grinding jar holder of the MM 301, the metallic parts can cause damage to the ceramic parts.		
Setting the frequency (vibration intensity)	If overloaded by a too high frequency setting, the display will alternate between the max. reached and the frequency set by you. Reduce the frequency setting to the max. value reached.			
Grinding process	Do not open the Plexiglas cover during the grinding process.	The analysis sets or grinding jars are brought to a standstill immediately by the built-in brake, but the grinding process cannot be continued for the remaining time. The MM 301 has to be restarted. The initial parameters are then once again available		
	Do not load the MM 301 unequally.	In order to ensure smooth running of the MM 301, both milling points must be loaded with the same mass.		
Cleaning	Disconnect the mains plug before every cleaning operation.	Danger of fatality from electric shock!		
	Do not clean the under running water	Danger of fatality from electric shock!		

ETSCH[®] International Distribution

tions in Milling & Sieving

Algeria

Retsch GmbH Rheinische Straße 36 42781 Haan Germany Phone: +49 (0)2129 5561-0 Fax: +49 (0)2129 8702 Email: info@retsch.com

Argentina

La Quimica Quirurgica S.A.C.I. Saavedra 247 C 1083 ACE Buenos Aires Argentina Phone: +54 115128 9100 Fax: +54 11 5128 9107 Email: abarie@laqq.com.ar

Australia

MEP Instruments PTY. Ltd. P. O. Box 1880 Unit 5, 112-118 Talavera Road NSW 2113 North Ryde Australia Phone: +61 2 9878 6900 Fax: +61 2 9878 6999 Email: info@mep-instruments.com.au Web: www.mep.net.au

Austria

Verder Ges.m.b.H. Eitnergasse 21 1230 Wien Austria Phone: +43 1 8651074-0 Fax: +43 1 8651076 Email: office@verder.at Web: www.retsch.at

Bangladesh

ENTECH LOGISTICS LTD. Plot#39, Block KA Pisciculture Society Shyamoli Dhaka - 1207 Bangladesh Phone: +880 2 8114789 Fax: +880 2 9122667 Email: entech@bdmail.net

Belarus

ODO "Belaquilon"

Laboratory Equipment Sharangovicha Str.19-217 220018 Minsk Belarus Phone: +375 172 590497

Fax: +375 172 590485 Email: info@belaquilon.by Belgium

Govatec NV

Kontichsesteenweg 17 2630 Aartselaar Belgium Phone: +32 3 8709640

Fax: +32 3 8709641 Email: info@retsch.be Web: www.retsch.be

Bolivia

Abendroth Intl SA POB: 3275 La Paz Bolivia Phone: + 591 2 418141 Fax: +591 2 2422855 Email: aisalab@entelnet.bo

Brasil

Altmann S. A. Blocko I - 70 andar 04794-000 Av. das Nacoes Unidas, 13.771 04794-000 Sao Paulo-SP Brasil Phone: +55 11 55073302 Fax: +55 11 55074196 Email: altmann@altmann.com.br

Brasil **IPAS Comercio e Representacoes**

C.G.C. 28.257.418/0001-06 Rua Conde de Lages No. 44 Sala 216 20241 Rio De Janeiro-Rj Brasil Phone: +55 21 2221 5893 Fax: +55 21 2259 4484 Email: ipasbr@attglobal.net Brasil

Polimate Ltda. Rua Basilio da Cunha, 206 Acclimacao 01544-000 Sao Paulo Brasil Phone: +55-11 5575 4041 Email: jose.neto@polimate.com.br

Bulgaria

MARVEL Ltd. 12, Kuklensko chaussee 4004 Plovdiv Bulgaria Phone: +359 32 678026 Fax: +359 32 678021 Email: sales@marvel.dir.bg

Canada Retsch Inc.

74 Walker Lane Newtown, PA 18940 United States Phone: +1 267 7570351 Fax: +1 267 7570358 Email: info@retsch-us.com Web: www.retsch-us.com

Chile

Equilab LTDA. El Quillay 627, sitio 96 Lampa Ciudad empr. Valle grande CASILLA 66 Santiago de Chile Santiago de Chile Chile Phone: +56 2 556 9974 Fax: +56 2 556 4006 Email: myalenzuela@equilab.cl

China

VERDER RETSCH (Shanghai) Trading Co.Ltd. Mira Garden Business Center Rm C 308 West Yanan Road 2633 200336 Shanghai China Phone: +86-21-61506046 Fax: +86-21-61023711 Email: info@retsch.cn Web: www.retsch.cn

Colombia

PC Quimica & Cia. Ltda Oficina 504, Edificio Agora Plaza Avenida 15 Nr. 124-31 Bogota, D.C. Colombia Phone: +571 612 5793

Fax: +571 612 5805 Email: cpadilla@multiphone.net.co

Colombia Symtek S.A.

(exclusivo para la industria de cementos y fundiciones) Cra 21 No 82-64/70 Of 702 Bogota Colombia Phone: +57 16167766 Fax: +57 16166889 Email: andres.diaz@symtek.com

Costa Rica

Tecnologia Aplicada Internacional POB: 1581-2150 San Juan Costa Rica Phone: +506 297 1011 Fax: +506 297 2231 Email: aulloa@tecnologia-aplicada.com

Croatia

Verder Ges.m.b.H. Eitnergasse 21 1230 Wien Austria Phone: +43 1 8651074-0 Fax: +43 1 8651076 Email: office@verder.at Web: www.retsch.at

Czech Republic

BIA VERDER s.r.o. Hlavni sklad Chlumecka 15 19800 Praha 9 - Kyje Czech Republic Phone: +420 286 853870 Fax: +420 286 856337 Email: jaroslav.seiml@retsch.cz Web: www.retsch.cz

Denmark

SKANLAB ApS Kvinderupvei 30 3550 Slangerup Denmark

Phone: +45 4738 1014 Fax: +45 4738 1064 Email: retsch@skanlab.com Web: www.skanlab.com

Ecuador

Quifatex SA Apartado 17-04-10455 10 de Agosto 10640 y Man. Zambrano Quito Ecuador Phone: +593 22477400 Fax: +593 4289734 Email: daniel.valladares@guifatex.com

Eavpt

Salzburg Co. for Scientific Equipment 292 Emtedad Ramsis 2

Nasr City - Cairo Egypt Phone: +202 22608082 Fax: +202 24041940 Email: salzburg_coeg@link.net

Estonia

G. W. Berg Oö Laki 25 12915 Tallinn Estonia Phone: +372-6505658 Fax: +372-6-505653 Email: laumets@gwb.ee Web: www.gwb.fi

Finland

OY G.W.Berg & Co. Ab MÑ kituvantie 7 01510 Vantaa Finland

Phone: +358 201 255255 Fax: +358 201 255256 Email: gwb@gwb.fi Web: www.gwb.fi

France VERDER S.A.R.L.

Parc des Bellevues Rue du gros chene 95610 Eragny sur Oise France Phone: +33 1 34643111 Fax: +33 1 34644450 Email: info@retsch-fr.com Web: www.retsch-fr.com

Germany

Retsch GmbH Rheinische Straße 36 42781 Haan Germany Phone: +49 (0)2129 5561-0 Fax: +49 (0)2129 8702 Email: info@retsch.com

E E E E C International Distribution

olutions in Milling & Sieving

Pakistan H. A. Shah & Sons

Office # 11, 3rd Fl., 109-W Sardar Begum Plaza Jinnah Avenue, Blue Area Islamabad Pakistan Phone: +92 51 2273557 Fax: +92 51 2272743 Email: hashah@apollo.net.pk

Peru

Quimica Suiza S.A. Apdo. 3919 2577 Av. Republica de Panama Lima Peru Phone: +51 1 2114055 Fax: +51 1 2114050 Email: atapia@quimicasuiza.com

Philippines

ITS Science & Medical Pte., Ltd.

219, Henderson Road, #011-02 Henderson Industrial Park Singapore 159556 Singapore Phone: +65 6273 0898 Fax: +65 6273 0810 Email: retsch@its-asia.com Web: www.lts-sciencemedical.com

Verder Polska Sp. z o.o.

Poland

ul, Ligonia 8/1 40-036 Katowice Poland Phone: +48 32 7815032 Fax: +48 32 7815034 Email: verder@verder.pl Web: www.verder.com

Portugal

Gravimeta Lda. Rua da Vilarinha, 1235 4100-517 Porto Portugal Phone: +351 226184232 Fax: +351 226184619 Email: geral@gravimeta.pt Web: www.gravimeta.pt

Republic of Korea (South)

Shinhan Science & Technology Samboo Bldg., 3rd Floor 537-6, Yatap-dong, Bundang-ku 463-840 Sungnam-si Kyunggi-do Republic of Korea (South) Phone: +82 31 706-3642 Fax: +82 31 706-3645 Email: scientech@kornet.net Web: www.shinhansci.com

Romania

Verder Romania S.R.L. Sector 5 Soseaua Villor no. 79 050152 Bucuresti Romania Phone: +40 21 3354592 Fax: +40 21 3373392 Email: eugen.ciotaru@verder.ro Web: www.verder.ro

Russia

RVS Co. Ltd. Khimanalit Cri. ul. Bumazhnaya, 17 190020 St. Petersburg Russia Phone: +7 812 7869516

Fax: +7 812 2520136 Email: retsch@rvs-ltd.ru Web: www.rvs-ltd.ru

Saudi-Arabía National Scientific Company Ltd Al Salam Building, Thrid Floor Tahlia Street Jeddah 21431 Saudi-Arabia Phone: +966 26644466 Fax: +966 26696192 Email: info@nsc-ksa.com Web: nsc-ksa.com

Singapore ITS Science & Medical Pte., Ltd.

219, Henderson Road, #011-02 Henderson Industrial Park Singapore 159556 Singapore Phone: +65 6273 0898 Fax: +65 6273 0810 Email: retsch@its-asia.com Web: www.its-sciencemedical.com

Slovak Republic

BIA VERDER s.r.o. Hlavni sklad Chlumecka 15 19800 Praha 9 - Kyje Czech Republic Phone: +420 286 853870 Fax: +420 286 856337 Email: jaroslav.selml@retsch.cz Web: www.retsch.cz

Slovenia

Verder Ges.m.b.H. Eitnergasse 21 1230 Wien Austria Phone: +43 1 8651074-0 Fax: +43 1 8651076 Email: office@verder.at Web: www.retsch.at

South Africa

Monitoring & Control Laboratories (Pty) Ltd. Ground Floor Lumley House 177 Jan Smuts Avenue Parktown North South Africa Phone: +27 11 327 6524 Fax: +27 11 327 6530 Email: gdavey@moncon.co.za Web: www.moncon.co.za

Spain

Biometa Tecnologia y Sistemas, S.A. Parque Tecnologico de Asturias Parcela 38

33428 Llanera (Asturias) Spain Phone: +34 902 244343 Fax: +34 985 269169 Email: info@biometa.es Web: www.biometa.es

Sweden

AB Ninolab Instrumentvaegen 4 Postfach/P.O. Box 137 19422 Upplands Vaesby Sweden Phone: +46 8 59096210 Fax: +46 8 59096210 Email: info@ninolab.se Web: www.ninolab.se

Switzerland

Schieritz & Hauenstein AG Finkelerweg 32 4144 Arlesheim Switzerland Phone: +41 61 7011060 Fax: +41 61 7011641 Email: schieritz@magnet.ch

Syria VART P.O.Box 4859 Damascus Syria Phone: +963 11 334 5858 Fax: +963 11 334 5857 Email: vart@mail.sy

Taiwan

Sunway Scientific Corporation 6F, No. 23, Lung Chiang Road 104 Taipei Taiwan Phone: +886 2 2771-8337 Fax: +886 2 2741-4646

Email: retsch@mail.sun-way.com.tw Web: www.sun-way.com.tw Thailand

Scientific Promotion Co., Ltd. 1759 Soi Wachirathamsathit 57 Sukhumvit 101/1 Bangchak Prakhanong, Bangkok 10260 Thailand

Phone: +66 2332 7960-4 Fax: +66 2331 8809 Email: marketing@spcgroup.co.th Web: www.spcpromotion.com

Tunisia

Retsch GmbH Rheinische Straße 36 42781 Haan Germany

Phone: +49 (0)2129 5561-0 Fax: +49 (0)2129 8702 Email: info@retsch.com

Turkey

KUTAY LABORATUAR CIHAZLARI TIC. A.S.

Levent Cad. Krizantem Sok. No:62 34330 1. Levent/Istanbul Turkey

Phone: +90 212 2829700 Fax: +90 212 2827975 Email: aysun@kutaygroup.com Web: www.kutaygroup.com

Ukraine

MACROLAB Ltd. Prospectus of 40-years of October 122/1 F. 91 03127 Kiev Ukraine Phone: +380 44 2583402 Fax: +380 44 2583402 Email: zhuravlov@voliacable.com Web: www.retsch.ru

United Arabian Emirates

Al Muhairi Scientific and Technical Supplies Est. 804,8th Floor, Green Tower, Deira P.O. Box 55219 Dubai United Arablan Emirates

Phone: +971 4 2623577 Fax: +971 4 2974560 Email: amsats07@eim.ae

United Kingdom

Retsch U.K. Ltd. Whitehouse Street, Hunslet Leeds LS10 1AD United Kingdom Phone: +44 113 2202666 Fax: +44 113 2465649 Email: info@retsch.co.uk

United States

Retsch Inc. 74 Walker Lane Newtown, PA 18940 United States Phone: +1 267 7570351 Fax: +1 267 7570358 Email: info@retsch-us.com Web: www.retsch-us.com

Uruguay

Alton S.A. Gral. Urquiza 3253 11600 Montevideo Uruguay Phone: +598 2 4872126 Fax: +598 2 487 6865 Email: altonventas@adinet.com.uy

page 3



Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

FAST SHIPPING AND DELIVERY

- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED
 SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment We also offer credit for buy-backs and trade-ins **www.artisantg.com/WeBuyEquipment** *>*

Instra View REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at **www.instraview.com** *>*

LOOKING FOR MORE INFORMATION?

Visit us on the web at <u>www.artisantg.com</u> for more information on price quotations, drivers, technical specifications, manuals, and documentation