14 423 00a

IKA®

IKA A 10

A10 072009

Reg.-No. 4343-01

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Operating instructions EN 7

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DE

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Warranty

In accordance with IKA warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

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Contents

EN

	Page
Warranty	6
Safety instructions	7
Correct use	7
Unpacking	8
Commissioning	8
Motor Circuit Breaker	8
Description	8
Grinding	8
Grinding while cooling	8
Replacing the grinding tool	9
Maintenance and cleaning	9
Accessories	10
Technical data	10
List of spare parts A 10	15
Spare parts diagram A 10	15

Safety instructions

In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.

Only staff who have been trained accordingly, know the appliance and are authorised to carry out work in this field should use this appliance.

The machine may only be opened by trained specialists - even during repairs. The machine is to be unplugged from the mains before opening. Live parts inside the machine may still be live for some time after unplugging from the mains.

The appliance can only be disconnected from the mains supply by pulling out the mains plug or the connector plug. The socket for the mains cord must be easily accessible.

Please follow the relevant safety instructions and guidelines, and occupational health and safety regulations for use in the laboratory. When working with the mill, the user must wear his personal protective equipment. Defective or inappropriate protective equipment can expose the user to the risk of spilling substances.

Please make certain that the mill is not tilted when it is in operation and hold it in place with both hands.

Please ensure that there are no parts of the body between the casing parts that move together.

Caution: Before removing the grinding beaker, switch off the mill and wait until the motor stops. Risk posed by trailing crushing elements (beaters or cutters).

Always switch off the appliance and disconnect the plug before cleaning or replacing parts.

Caution: The grinding beaker, crushing elements and rotor shaft may be hot after long grinding cycles. Allow the mill to cool down before cleaning or before changing crushing elements.

Caution: Beware of the sharp - edged cutter in particular when changing crushing elements and when cleaning. Always wear

your personal protective equipment (protective gloves, etc.). Please note that wear debris from the appliance or rotating accessories can end up in the medium to be processed. Please contact IKA should you have any questions in this respect.

The analytical mill **A 10** may not be used in explosive atmospheres or under water. Do not crush explosive, toxic or noxious substances. Please observe the relevant safety instructions and guidelines regarding dust explosions

Always wear suitable protective equipment when grinding the mill feed with coolants (liquid nitrogen or dry ice), (see laboratory guidelines and accident prevention regulations applicable to the laboratory).

Ensure that pressure does not build up in the grinding chamber (particularly when using coolants).

Caution: When using liquid nitrogen, you must only attach on the grinding beaker and begin the grinding process when the liquid nitrogen is completely vaporized.

Caution: Appliance parts may still be extremely cold for some time after cooling with liquid nitrogen and should only be handled whilst wearing protective equipment.

Correct use

IKA analysis grinders A 10 is suitable for carrying out low-loss grinding of hard and brittle substances as of grain size 6 - 7mm. The grinders operate without creating dust or lost. Beaters, blades and grinding are replaceable. The useful volume of the grinding chamber is 50 ml for the A 10 grinder.

Hard and brittle ground material is crushed by means of beaters or blades rotating at high speed.

Oil substances or those that are given to smearing and which are not capable of being ground at room temperature have to be cooled for grinding (e.g. brittle with liquid nitrogen).

Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder). The A 10 delivery quantity includes: one A 10 grinder, one A 14 stainless steel beater, one A 18 grinding chamber reduction element, one set of tools and operating instructions.

Description

The device consists of a drive unit with grinding chamber an cover. They are connected to each other by means of quick-release locks. When assembling the two units make certain that the safety pin in the drive unit is pushed down by the pressure pin on the cover. Otherwise, a tight and reliable connection is not ensured, the yellow error lamp will light up and the grinder will not start. The device is started by pressing the push-button at the front.

Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged.

Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.

Note:

The grinder is designed for short term operation. Maximum operating time: A 10 - 5 minutes .

After that, the grinder must cool off for approximately 10 minutes. This condition does not apply if the grinders is connected to a cooling unit.

Motor Circuit Breaker

The device is electronically protected against overload and excessive temperatures. A error yellow lamp lights on the front panel up if the fuse is triggered and the motor comes to a standstill.

Grinding

Fineness is mainly dependent on how long the grinding operation is performed.

Dry, brittle materials are easily ground, but may not be larger than the size of a pea as there may be danger of damage when carbide metal beaters are used. The A 14 stainless steel beater is assembled from rust- and acid-proof stell (material no. 1.4034). The A 15 carbide metal beater made of carbide metal (86% tungsten carbide and 14% cobalt) is available if particulary wear-resistant material is required. The A 17 star-type blade (material no. 1.4310) is available for specifically light ground materials such as pharmaceuticals, tablets, cardboard, felt, polystyrene, plastic foams (Brittle) and the like

The A 18 grinding chamber reduction element (included in the delivery quantity) is inserted into the grinding chamber for comminution of small sample quantities (1 - 5ml) and specifically lighter products.

Grinding while cooling

Plastic, soft or smeary materials have to be cooled if they are to be ground.

There are two hose couplings located at the rear of the grinder in

order to connect the cooling brine. Cooling with tap water is usually sufficient for substances which become hot as a result of the grinding operation and are thus given to smearing or those that produce oil.

The grinding chamber may only be filled halfway especially when drugs containing oils and similar materials are ground.

If more refrigeration is required then the ground matter should be cooled prior to being put into the grinding chamber.

The following procedure shall apply if liquid nitrogen is used in order to embrittle the material to be ground:

Safety infomation:

Suitable protective gloves must be worn when filling and transferring liquid nitrogen!

Protective glasses must be worn!

- The sample to be ground is put into a suitable vessel (e.g. evaporating dish) and covered with enough liquidnitrogen until it is completely covered.

When using liquid nitrogen make certain that the sample is not capable of reacting with oxygen.

- Cool grinding chamber using liquid nitrogen.
- Put the deep-cooled sample into the grinder and grind approximately 15 seconds.
- Then open the grinder and add liquid nitrogen for only as long as it takes for the sample to become deep-cooled again (embrittled).

Safety infomation:

Make certain that the liquid nitrogen has evaporated before beginning the grinding operation.

Note: Danger of overpressure!

- Repeat the grinding operation for approximately 15 seconds.
- Repeat the sequence of cooling and grinding until desired fineness has been achieved.
- If the required fineness has not been achieved after 5 minutes, then a further reduction in grain size is unlikely.

The respective grinding chambers are sealed by means of labyrinth sealings. Though they are dust-proof, they are not leak-proof.

This should be taken into account for substances where liquids may escape.

Replacing the grinding tools

The grinding tool is fixed by means of the retaining tube and the cap nut is unscrewed by using zhe SW 8 or 10 spanner. The grinding tool, which rests on a double-edge track, may then be pulled out.

Maintenance and cleaning

The IKA A 10 is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate. When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices.

To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices.

Electrical equipment may not be put into the cleansing agent for cleaning.

The operator shall consult the manufacturer before using a method of cleaning or decontamination other than the one which is re-

commended by the manufacturer in order to make certain that the method which the operator intends to employ will not damage the equipment.

Accessories

- A 18** grinding chamber reduction element
- A 11** cover
- A 14** spare beater, stainless steel
- A 15** beater, carbide metal
- A 17** star-type blade, stainless steel

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
List of spare parts

Item	Designation	Item	Designation
7	cup nut	29	threaded bushing
9	grooved pin	31	holder
13	distance ring	32	foundation
14	O-ring	34	font plate
15	O-ring	39	spare cutter A14
16	Nilos-ring	40	ginding bowl
17	Nilos-ring	43	spring
19	ball bearing compensating disc	45	eccentric lever
20	fan blade	48	commutator motor
21	cap	50	printed circuit board
22	protective rubber sleeve	51	cover A11
25	housing	52	connection line
26	flange	2001	grooved ball bearing
27	felt ring	2002	carbon brush
28	hub	3001	push-button
		3002	push-button

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China

IKA® Works Guangzhou
173 - 175 Friendship Road
Guangzhou
Economic and Technological
Development District
510730 Guangzhou, China
Tel. +86 20 8222-6771
Fax +86 20 8222-6776
sales@ikagz.com.cn

Asia
Australia

IKA® Works (Asia)
Sdn Bhd
No. 17 & 19, Jalan PJU 3/50
Sunway Damansara
Technology Park
47810 Petaling Jaya
Selangor, Malaysia
Tel. +60 3 7804-3322
Fax +60 3 7804-8940
sales@ika.com.my

Japan

IKA® Japan K.K.
293-1 Kobayashi-cho
Yamato Koriyama Shi, Nara
639-1026 Japan
Tel. +81 743 58-4611
Fax +81 743 58-4612
info@ika.ne.jp

Korea

IKA® Korea Co LTD
1710 Anyang Trade Center
1107 Buhung-dong,
Dongan-gu
Anyang City, Kyeonggi-do
Post code: 431-817
South Korea
Tel. +82 31-380-6877
Fax +82 31-380-6878
michael@ikakorea.co.kr

India

IKA® India Private Limited
814/475, Survey No. 129/1
Mysore Road, Kengeri
Bangalore - 560 060
Karnataka, India
Tel. +91 80 26253900
Fax +91 80 26253901
info@ika.in

Brasilia

IKA® Works Inc.
Rua São Bento, 701 Sala 1
CEP 13160-000 Centro -
Artur Nogueira - SP
Brasil, South America
Tel. +55 19 3877 2399
Fax +55 19 3877 2399
fcabral@ika.net