IKA RW 20.n
IKA RW 20 DZM.n
IKA RW 20 DZM.n P4

BETRIEBSANLEITUNG	D	4
OPERATING INSTRUCTIONS	GB	10
MODE D'EMPLOI	F	16
INSTRUCCIONES DE MANEJO	E	22
HANDLEIDING	NL	28
INSTRUZIONI PER L'USO	I	34
DRIFTSANVISNING	S	40
DRIFTSINSTRUKS	DK	46
KÄYTTÖOHJE	N	52
DRIFTSVEJLEDNING	SF	58
BETRIEBSANLEITUNG	GR	64
INSTRUÇÕES DE SERVIÇO	P	70

RW201095EU

Reg. No. 4343-01


Fig. 1


Fig. 2


Fig. 3

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D

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Contents

	Page
Safety instructions	10
Correct use	11
Unpacking	11
Useful facts	11
Drive	11
Gear speed change	12
Motor circuit breaker	12
Speed display	12
Commissioning	13
Switching on the machine	13
Output shaft	13
Maintenance and cleaning	13
Technical data	14
Accessories	15
Permitted IKA stirrer tools	15
Associated standards and regulations	15
Guarantee	15
Spare parts list	82
Spare parts diagram	83/84

Safety instructions

Never operate the equipment with the stirrer tools rotating freely. Ensure that parts of the body, hair or items of clothing cannot be trapped by the rotating parts. Pay attention when setting the speed to any imbalance of the stirrer tools and possible spraying of the medium to be stirred. Use a stirrer shaft protection device!

The equipment is not suitable for manual operation. Please follow the relevant safety instructions and guidelines, and occupational health and safety regulations for use in the laboratory.

Avoid knocks and impacts on the lower end of the shaft and the chuck gear teeth. Even minor, invisible damage can lead to imbalance and uneven shaft action. Careful handling guarantees safe working and the long life of the equipment.

Imbalance of the output shaft, the chuck and in particular the stirring tools can lead to uncontrolled resonant vibrational behaviour of the equipment and the whole assembly. Glass apparatus and stirrer containers can be damaged or shattered by this. It can cause injury to the operator, as can the rotating stirring tool. If vibration of the equipment is noticed, the speed must be reduced immediately in all cases until no more vibration occurs. In this case exchange the stirring tool for one without imbalance or remedy the cause of the imbalance.

Operating with a freely rotating shaft end is dangerous. Therefore for safety reasons the stirrer tool is permitted to project over the upper edge of the housing only when the machine is not running.

The unit must be switched off before changing to another gear speed, otherwise the toothed gears in the reduction stage will be damaged.

The machine may only be opened by trained specialists - even during repairs. The machine is to be unplugged from the mains before opening. Live parts inside the machine may still be live for some time after unplugging from the mains.

The high torque developed by the RW 20 DZM.n P4 requires particular care in the choice of stand, cross sleeve and anti-rotation element for the stirring bowl.

Correct use

Stirrers in the RW20.n series are suitable for stirring and mixing liquids with low or high viscosity. They are designed for use in the laboratory. For correct use the machine must be fixed to a stand.

possible to work on stands - even with high glass apparatus under vapour extraction hoods - without extension leads. The precision stirrer shaft is designed as a hollow shaft to allow the stirrer shanks to be inserted.

Drive

Unpacking

Please unpack the machine carefully and inspect for damage. It is important that any transit damage should be noted at the time of unpacking. In certain circumstances it may be necessary to investigate immediately (post, rail or freight forwarder). The guarantee covers: One RW20(.) stirring machine, one extension arm, one hexagonal socket screw, one hexagon socket offset screw key and operating instructions.

The method of operation of the motor is characteristic of the operating behavior of the RW 20.n series units. The motor can be infinitely adjusted via a friction wheel drive unit but is operated in one position only. The power output, rotational speed and torque of the motor are regarded as constant and are optimized for this operating position. After the friction wheel drive the motor power output is transmitted to the off-drive shaft of the stirrer via either the first or second stage of the subsequent manually adjustable toothed gear train.

Neglecting all losses, the power at the off-drive shaft during power transmission (conversion) is always as great as the power output of the motor. The gear train merely provides rpm and torque conversion, which can be described in purely theoretical terms by the characteristic curve plotted in Fig. 4.

Useful facts

In buying this machine you have chosen a high quality product. Ease of use is guaranteed with its handy shape and simple operation. The housing offers protection against ingress of liquids. The materials used and their precise identification make recycling possible and simple, and enables re-use of the parts. The speed can be infinitely adjusted between 60 and 2000 rpm 50 Hz; 72 and 2400 rpm 60Hz. The machine can be operated continuously and warms up during operation. The generously proportioned cooling surfaces enable distribution and transference of heat to be as even as possible. The two metre long mains lead makes it

Friction wheel wear is kept low by a helical clutch matched to the gear train. The down-pressure required at the friction wheel is adjusted by the helical clutch according to the moment acting on the stirrer shaft. A low torque results in a low down-pressure; a high torque results in a high down-pressure.

In order to be able to process high-viscosity media, the RW 20 DZM.n P4 does not have the second switchable reduction stage, but features instead a planetary gear train with a gear ratio of 3.7 : 1 immediately after the first toothed gear stage. This provides a speed range of 16 to 135 rpm at 50 Hz and 20 to 160 rpm at 60 Hz.

Gear Speed Change

RW 20.n and RW 20 DZM.n only

(See Fig. 5)

Should it be necessary to change the gear speed in the course of setting up an experiment or process with glass apparatus, then attention must be paid to the movement of the stirrer attachment by the eccentric bearing bush and to the movement in the axial direction. If this is unacceptable, the stirrer attachment must be loosened in the chuck.

Turning (counter-clockwise) the bearing bush at the black knurled surface over the chuck disengages the gears. The center axis of the stirrer attachment shifts 2.2mm in the radial direction. The bearing bush is then moved 5.5mm up or down in the axial direction according to the range to be set. Turning back (clockwise) the bearing bush re-engages the gears and swings the center axis of the stirrer attachment back to its original position.

Motor Circuit Breaker

The ball race capacitor motor, protected to IP 42, is maintenance-free. In its winding it carries a special overheating circuit breaker (self-locking temperature limiter) which shuts down the motor permanently if the permissible motor temperature is exceeded.

This special type of circuit breaker is used when the automatic reactivation of the protected unit after overheating is undesirable or impermissible. The principle of the temperature limiter is that of a switch that heats up and isolates the motor from the power supply when the maximum permissible temperature is reached.

When the self-locking temperature sensor is tripped, the yellow glow lamp (**B**) on the front plate of the unit lights up to indicate the motor overload and the stoppage of the motor. In this case the unit can only be switched back on after it has been switched off for approx. 4-5 minutes by pressing the rocker switch (**A**) on the front plate or by disconnecting the mains plug.

After this short cooling period, the contacts in the current circuit are closed again and the unit is operational once more.

Speed Display

The rotation speed is adjusted with the rotary knob (**C**) on the front plate of all RW 20.n stirrer units.

RW 20.n

In the RW 20.n unit, a defined rotational speed is set via the scaling of the rotary knob position. The plate (**Fig. 5**) attached to the side of the unit lists the speeds attained at the set scale values. (The values apply only when low loads are acting on the stirrer shaft.)

RW 20 DZM.n

In the RW 20 DZM.n unit, the rotational speed is displayed directly in revolutions per minute (rpm) on the LCD display (**D**). The display can only show speeds up to max. 1999 rpm, which means that the off-drive shaft may be moving at a faster speed than displayed when there are low loads acting on the stirrer shaft.

RW 20 DZM.n P4

In the RW 20 DZM.n P4 unit, the rotational speed is displayed directly by the LCD display (**D**), as in the RW 20 DZM.n. As the toothed gear train is no longer adjustable and is set to the lowest speed, the entire rpm range can be displayed.

Commissioning

For correct use, the stirrer machine must be fixed with a cross sleeve (eg R181) to a stable stand (eg R1822). The stirrer container should always be fixed with a tensioner (eg RH1) for safety reasons.

Switching on the machine

Check whether the voltage given on the type plate corresponds to the available mains voltage. The socket used must be earthed (fitted with earth contact). If these conditions have been met, the machine is ready to operate when the mains plug is plugged in. If these conditions are not met, safe operation is not guaranteed and the machine could be damaged.

During commissioning of the machine the output shaft starts to run at the last speed set. Therefore check the setting of the control knob. Also ensure that the speed set is suitable for the test texture selected. If in doubt, set the speed knob (B) to the lowest speed (left-hand position). In case of doubt, turn the speed knob (C) to the lowest speed (all the way left) and set the gear train to the lowest reduction stage.

The unit starts to operate when the rocker switch is pressed.

Note: **The rocker switch (A) is not a mains power switch !**

Output shaft

If the unit has not been used for some time, a knocking noise will be heard when switching on, which is due to the preload on the friction lining of the friction wheel drive. This has no detrimental

effect on the function of the unit, and the knocking noise will disappear after a short running-in period.

The chuck and off-drive shaft allow stirrer attachments with shank diameters of up to 10mm to be clamped and inserted (except RW20DZM.n P4). The opening at the top of the housing allows stirrer attachments to be pushed up above the top edge of the housing, e.g. in order to change the stirring bowl with the unit switched off.

Maintenance and cleaning

The stirrer is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Only use water with a detergent additive containing a surfactant for cleaning, or for heavier soiling isopropylalcohol.

When cleaning the unit, take particular care not to bring substances containing oil or grease into contact with the running surfaces of the cone pulley and the friction wheel. This would substantially reduce the frictional coefficient resulting from the pairing of the friction wheel and cone pulley materials, and would impair power transmission in the unit.

Technical data

Speed range (RW20.n / RW20 DZM.n)

(50Hz stage I)	min⁻¹	60 - 500
(50Hz stage II)	min⁻¹	240 - 2000
(60Hz stage I)	min⁻¹	72 - 600
(60Hz stage II)	min⁻¹	288 - 2400
(RW20 DZM.n P4)		
(50Hz stage I)	min⁻¹	16 - 135
(50Hz stage II)	min⁻¹	20 - 160

Max torque stirrer shaft

(RW20.n / RW20 DZM.n)	Ncm	185
(RW20 DZM.n P4)	Ncm	700

Permitted on-time:

Toggle knob

Speed display: **(RW20.n)** Scale on front plate of unit
(RW20 DZM.n and RW20 DZM.n P4) LCD display

Measurement fault: max. ±0,5% ±30 Digit

Nominal voltage: **VAC** 230±10%

VAC 115±10%

Frequency: **Hz** 50 EURO

Hz 60 USA

Input power: **W** 72 EURO

W 87 USA

Power output: (short term)
 (constant operation) **W** 35

W 20 +35%

Overall efficiency: **%** 40

IP 20

Protection type to DIN 40 050:
 Exess voltage category: **II**

Contamination level: **2**

Protection at overloaded:
 Temperature sensor in motor winding

A 4 T (IKA ident.no.

25 851 00)

Ambient temperature: **°C** +5 to +40

% 80

Ambient humidity: (rel.)

Operating position:

on stand, clamping chuck pointing down

Drive:

Rib-cooled capacitor motor with
 friction wheel drive and subsequent
 2-stage toothed gear train.
 (RW 20 DZM.n P4 fixed at stage I;
 planetary gear train i=3.7:1)

Maximum stirring quantity - water:

ltr 20

For viscosity: **(RW20.n / RW20 DZM.n)** MV
(RW20 DZM.n P4) HV

Viscosity	mPa · s	
VL very low viscosity	0 100	Water up to thin-bodied oil
LV low viscosity	100..... 1000	Thin-bodied oil up to thick-bodied oil
MV medium viscosity	1000.. 10000	Thick-bodied oil up to honey at approx. 20 °C
HV high viscosity	10000 100000	Honey at approx. 20 °C up to inks/dyes

Clamping chuck clamping range: **mm** 0,5 -10

Hollow shaft internal diameter: **mm** 10,5

Etension arm:

(RW20.n / RW20 DZM.n) **mm** 13x160 long

(RW20 DZM.n P4) **mm** 16x160 long

Housing: Al diecasting, thermoplastic

	RW20.n	RW20 DZM.n	RW20 DZM.n P4
Dimensions without exten- sion arm (WxDxH) mm	88x188x292	88x210x292	88x210x402
Gewicht mit Ausleger und Spannfutter: kg	2,9	3,1	4,6

Accessories

RW 20.n and RW 20 DZM.n

R 1823 Plate stand

R 181 Cross sleeve

R 182 Cross sleeve

RW 20 DZM.n P4

R 2723 Telescopic stand

R 201 Cross sleeve

R 271 Cross sleeve

RH 1 Tensioner

FK 1 Flexible coupling

RH 2 Tensioner

R 263 Bowl clamp

RW 20 DZM.n P4

R 1343 Propeller stirrer 4fl

1 000

R 1345 Propeller stirrer 4fl

800

R 1331 Anchor stirrer

1 000

R 1393 Worm stirrer

1 000

Permitted IKA stirrer tools

Max. Speed 1/min

For all options

R 1313 Turbine stirrer 800

R 1375 Surface stirrer 800

R 1330 Anchor stirrer 1 000

RW 20.n and RW 20 DZM.n

R 1342 Propeller stirrer 4fl 2 000

R 1381 Propeller stirrer 3fl 2 000

R 1382 Propeller stirrer 3fl 2 000

R 1389 Propeller stirrer 3fl PTFE 800

R 1311 Turbine stirrer 2 000

R 1312 Turbine stirrer 2 000

R 1300 Dissolver stirrer 2 000

R 1303 Dissolver stirrer 2 000

R 1352 Centrifugal stirrer 2 000

R 1373 Surface stirrer 2 000

R 1332 Anchor stirrer PTFE 800

R 1335 Kneading stirrer 2 000

Associated standards and regulations

Associated EU guidelines

EMV guidelines: 89/336/EC

Machine guidelines: 89/392/EC

Construction in accordance with the following safety standards:

EN 61 010-1 / VDE 411-1 EN 60 204-1

EN 50 081 EN 292-1, -2

EN 50 082-1 EN 414

EN 55 011 IEC 1010-1

EN 55 014-1 UL 3101-1

EN 60 555-2, -3 CAN/CSA C22,2 (1010-1)

Guarantee

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 12 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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List of spare parts






• RW 20.n + RW 20 DZM.n .. RW 20 DZM.n + RW 20 DZM.n P4 ... RW 20 DZM.n P4

Item	Designation	Item	Designation
4	front plate	62 ..	intermediate shaft
5	housing	63 ..	distance tube
6	motor support	63 ...	threaded pin
7	A/C-motor	64 ..	grooved ball bearing
11	insulating hose	65 ..	eccentric bearing
13	flat spring	66 ..	off-drive shaft
14	conveyor plate	67 ..	perforated washer
15	adjusting spindle	68 ..	pressure spring
18 ..	plate	68 ...	grooved ball bearing
22 ..	insulating plate	69 ..	steel ball
23	cone pulley	69 ...	grooved ball bearing
24	fixing angle	70 ..	threaded pin
25	Sechskantmutter	71 ..	hexagonal nut
26	glow lamp	71 ...	gear-wheel rim chuck
27	switch	72 ..	cover plate
29	toggle knob	73 ..	retaining ring
31	insulating plate	74 ..	gear-wheel rim chuck
32	screw	75 ..	threaded pin
33	protective rubber sleeve	75 ...	planetary gear train
34	cable clamp	76 ..	ball bearing shim
35	PCB mains fuse	76 ...	distance tube
37	flange	77 ..	perforated washer
38	friction wheel	78 ..	threaded pin
39	clutch member	79 ..	flange
40	pressure spring	80 ..	bearing flange
41	protection against torsion	80 ...	bearing flange
51 ..	housing	81 ..	bearing housing
52 ..	distanzrohr	82 ..	distance tube
53 ..	window	83 ..	packing
54 ..	shaft	86 ..	intermediate shaft with gear
55 ..	PCB display	87 ..	toothed gear shaft
56 ..	PCB mains	2001	gear unit (replacement)
58 ..	measuring sensor	3001 ..	grooved ball bearing
59 ..	connecting cable	5001	grooved ball bearing
61 ..	off-drive shaft	5002	fan rotor
		5003	fan cowl

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


IKA® WORKS, INC.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

The Global Market of Laboratory Equipment
<http://www.labworld-online.com>

Europa - Afrika

IKA® -WERKE GMBH & CO.KG
Janke & Kunkel-Str. 10
D 79 219 STAUFEN
GERMANY
TEL. 07633/831-0
FAX 07633/831-98
E-mail: sales@ika.de
Internet: <http://www.ika.net>

America

IKA® WORKS, INC.
2635 NORTH CHASE PKWY. SE
WILMINGTON, NC 28405-7419
TEL. 800/733-3037
TEL. 910/452-7059
FAX 910/452-7693
E-mail: usa@ika.net

Asien - Australien

IKA® Works (Asia) Sdn Bhd
(Company No. 340448-K)
Lot 2, Jalan Indah 1/2
Taman Industri Rawang Indah
48000 Rawang
Selangor, Malaysia
TEL : (603) 6093 3322
FAX : (603) 6093 3940
E-mail: ika@tm.net.my