## **BIOVENDIS**



# Bead Ruptor 4

MINI BEAD MILL HOMOGENIZER

Powerful performance in a compact design

Built upon the innovative design and unsurpassed quality of the Omni Bead Ruptor family.

## Bead Ruptor 4

MINI BEAD MILL HOMOGENIZER

**ORDER NO. 25-010** 



#### **Processing Range:**

 $25 \,\mu L$  to  $6 \,mL$ 

#### **Power Rating:**

100-230 VAC/50-60 Hz

#### **Dimensions:**

8.5" W x 11.5" D x 10" H

Weight: 16.5 lbs. (7.48 kg)

Run Time: 1 - 5 mins.

Programs: 1

Cycles: 1

Performance Range:

 $30 \times g - 150 \times g$ 

#### **Efficiently lyse:**

- > Human tissues
- > Animal tissues
- > Plant & food
- > Micro-organisms



For best results, use with Omni pre-filled bead kits

### Designed for personal use

- > Broad performance range gentle mixing to high force homogenization
- > Small footprint requires minimal bench space
- > Holds 4 x 0.5 mL / 1.5 mL / 2 mL or 1 x 7 mL tubes

Patent Pending Snap In Tube Carriage



## Omni Bead Ruptor 4 is ideally suited for:

- > Homogenizing up to 4 samples in seconds
- Reproducible extraction of proteins and nucleic acids
- Simultaneously processing 4 x 0.5 mL, 4 x 1.5 mL, 4 x 2 mL or 1 x 7 mL samples

### > For more information — www.biovendis-products.de

#### Biovendis Ltd.

Am Exerzierplatz 3; D-68167 Mannheim **EMAIL** info@biovendis.de **TEL** 0621 - 43 00 52 52 **FAX** 0621 - 43 00 52 50

## **BIOVENDIS**

Biovendis is a partner of Biolabproducts GmbH

