
KERN **572-37**







Measuring system	
Adjustment options:	External calibration
Eccentric loading at 1/3 [Max]:	50 mg
Linearity:	50 mg
Readability [d]:	10 mg
Recommended adjusting weight:	2 kg (F1)
Repeatability:	20 mg
Resolution:	300000
Stabilisation time under laboratory conditions:	3 s
Warm up time:	2 h
Weighing capacity [Max]:	3 kg
Weighing system:	Strain gauge
Weighing units:	kg to tl (HK) g mo tl (Cn) lb tl (Singap. Malays) dwt ct tl (Tw) oz gn t

Display	
Display digit height (large):	1,800 cm
Counting	
Counting resolution:	300000
Mininum piece weight at piece counting (Laboratory:	10 mg
Mininum piece weight at piece counting (Normal):	100 mg

Functions	
Conversion function:	yes
Counting function:	yes
Percentage determination:	yes
PRE-TARE function:	yes
Vibration-free weighing (Animal weighing program):	yes

Maximum humidity:	80 %	
Maximum operating temperature:	40 °C	
Minimum ambient temperature:	10 °C	

Power supply	
Charging time:	10 h
Input voltage:	100 V - 240 V
Mains adapter:	Mains adapter external
Mains adapter/adapter included:	EURO CH US UK
Operating time (Backlight off):	30 h
Operating time (Backlight on):	15 h
Rechargeable battery:	Rchrg, battery optional

DAkkS Certificate:	963-127
Category	
Brand:	KERN
Category:	Balances
Product Group:	Precision balance
Packaging & shipping	

Services (optional)

r dokaging a shipping	
Delivery:	24 h
Dimensions packaging (WxDxH):	400 x 340 x 220 mm
Gross weight:	4,400 kg
Net weight:	2,300 kg
Shipping method:	Parcel service

Construction	
Casing material:	Cast aluminium alloys
Dimensions housing (WxDxH):	180 x 310 x 85 mm
Level indicator:	yes
Material weighing plate:	stainless steel
Overall dimensions mounted (WxDxH):	180 x 310 x 90 mm
Weighing surface (d):	150

KERN **572-37**



All-rounder e.g. as precision balance in the laboratory or in harsh industrial applications

Piktogramme

STANDARD



OPTION

